POLYFIX®

Bedienungsanleitung
Operating manual
Mode d'emploi
Istruzioni per l'uso
Instrucciones de manejo

Artikelübersicht - Product list - Liste de produit**Lista dei prodotti - Relación de productos**

Dispenser - Dispenser - Distributeur - Dosatore - Dispensador
POLYFIX®

101.070-38 2 - 10 ml
101.070-48 10 - 50 ml

Ventilblockeinheit - Valve block unit - Bloc soupape - Blocco valvole - Bloque de válvulas
POLYFIX®

101.075-38 2 - 10 ml
101.075-48 10 - 50 ml

Ersatzkolben - Spare piston - Piston de rechange - Stantuffo di ricambio - Émbolo de recambio
POLYFIX®

101.077-38 2 - 10 ml
101.077-48 10 - 50 ml

Verschlußkappe - Stopper cap - Capuchon de fermeture - Tapetto - Caperuza de cierre
rot - red - rouge - rosso - roja

101.079 10 ml + 50 ml

Volumen-Einstellsystem - Volume setting system - Réglage du volume - Sistema di regolazione del volume
- Sistema regulador del volumen

101.000-KS 10 ml + 50 ml

Adapter - Adapter - Adaptateur - Adattatore - Adaptador
PP

101.090-25 GL 25
101.090-28 GL 28
101.090-30 GL 30
101.090-33 GL 33
101.090-36 GL 36
101.090-38 GL 38
101.090-40 GL 40
101.090-45 GL 45
101.090-52 GL 52
101.090-55 GL 55
101.090-60 GL 60
101.090-65 GL 65

Kegelschliff-Adapter - Adapter ground joint - Adaptateur pour rodages conique - Adattatore smerigliato conico
- Adaptador esmerilado cónico

PP

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GB

Preface

We thank you for the purchase of a POLYFIX® bottle top dispenser. You did a good choice. With more than 20 years of experience in making dispensers we guarantee a high standard of product quality and safety.

To guarantee quality and safety of the dispenser also for a long term, a correct operating and maintenance is necessary.

Therefore please read this operating manual carefully before using the dispenser and take notice of the warning instructions as well as the corresponding regulations.

The delivery contents: 1 dispenser, 1 intake tube, 2 adapters made of PP, 1 operating manual. Please check this before using the dispenser.

Operating Manual Bottle Top Dispenser POLYFIX®

I. Safety Instructions

1. Before using the dispenser please read the operating manual carefully and follow the safety instructions.
2. Observe general notes on hazards and safety regulations, e. g. wear protective clothing, gloves and glasses.
3. Observe any hazard symbols on reagent bottles.

4. Use only appropriate containers to dispense liquids. If using small bottles, please use a stand.
5. Always make sure that the discharge tube does not point at the user or any other person.
6. When transporting the dispenser please ensure that the piston is in zero-position and no liquid is in the cylinder.

**The stopper cap must close the discharge tube and the volume setting system must be set at zero-position!
Never carry the dispenser at the finger guard, only at the PP valve block!**

7. Please check the correct assembling of the dispenser and the tight fit on the bottle before use.
8. To guarantee a trouble-free working please clean and use the dispenser according to the operating manual.
9. Repairs if necessary must be carried out by the manufacturer himself or authorized service stations only by using original parts.

II. The World of POLYFIX®

POLYFIX® 10 ml

Art. No. 101.070-38

Bottle top dispenser for aqueous liquids and acids. Piston made of precision ground glass. Integrated glass dispensing cylinder, plastic coated. Valve block made of PP. Valve system with Hastelloy metal spring in discharge channel. FIX adapter.

POLYFIX® 50 ml

Art. No. 101.070-48

Bottle top dispenser for aqueous liquids and acids. Piston made of precision ground glass. Integrated glass dispensing cylinder, plastic coated. Valve block made of PP. Valve system with Hastelloy metal spring in discharge channel. FIX adapter.

Warning:

Please see the technical information about polypropylene.

This dispenser must not be used with solvents!

III. Modern Engineering

Finger guard

Finger guard made of PP. All descriptions on the finger guard and the graduation on the scale are lasered to be indestructable.

Volume setting system

Easy-to-use volume setting system with fine adjustment for calibration. Suitable to be used by left- and right-handers.

Dispensing cylinder

Integrated glass cylinder, vacuum-calibrated and plastic coated, with stopper ring. Graduation amber stain diffused.

Valve block

Valve block made of PP. FIX Adapter.

Valves


Intake and discharge valve made of PP. Cristal valve balls. Hastelloy metal spring in discharge channel.

Dispensing piston

Dispensing piston made of precision ground glass.

Discharge tube

Integrated discharge tube made of PP. For safety reasons the discharge tube can not be taken apart. Red stopper cap made of PP.


IV. Handling

1. Assembling / Filling

Before using the dispenser please check the dispenser for any damages or glass breakage. Push intake tube over the hose connection of the PTFE valve block.


Screw dispenser direct onto containers with thread GL 32 or use adapters to fix on any other threads or ground joint necks. Check the correct fitting. The dispenser can be turned on the flask and fixed in any direction by using the built-in FIX adapter.

For filling the dispenser the first time lift piston several times just a few until no air bubbles are seen in the cylinder.

2. Volume setting

Loosen the nut of the volume setting system, set the desired volume. The upper edge of the fine adjustment screw shows the volume set. Fix the screw.

For calibrating purposes turn the screw of the fine adjustment. Turn downwards to increase the volume, turn upwards to decrease the volume.

Warning:
If the stopper cap is not removed when dispensing, there may be an over-pressure in the channel system and uncontrolled dispensing may be the consequence.

3. Dispensing procedure

- a) For safety make sure that the discharge tube does not point at the user or any other person. By using the built-in FIX adapter it is possible to fit the dispenser accordingly. Remove stopper cap by holding the discharge tube tight with other hand.
- b) Set the desired volume with the volume setting system.
- c) Lift the dispensing piston up with the finger guard. Do not touch the piston itself to avoid contact with chemicals. The piston stops at the volume set. The graduation on the glass cylinder is an additional optical control of the volume set or air-bubbles.
- d) To dispense the liquid push the piston down.
- e) Avoid uneven piston movements or hard pushes against the stopper ring.**
- f) Replace stopper cap after the dispensing procedure by holding the discharge tube tight with other hand.

Warning:

The piston must be pushed down manually. Therefore chemicals may remain in the dispensing cylinder after dispensing. To avoid unintentional dispensing, always dispense the total volume set. The stopper cap must be placed on the discharge tube after every use.

When dispensing observe safety precautions!

V. Operating Temperature

The operating temperature is from +15° C (+59° F) to +50° C (+122° F).

VI. Sterilization

The dispenser can be autoclaved at +121° C (+248° F) / 1 bar. The piston can remain in the cylinder when autoclaving.

Before sterilization loosen screw of the volume setting system.

Warning: Do not autoclave with higher temperatures. This will damage the dispenser. Wait before use until the dispenser has room temperature.

VII. Transport and Storage

- 1. When transporting the dispenser never carry on the finger guard but always on the PP valve block to avoid unintentional dispensing of liquid! Replace the stopper cap before transporting and fix the volume setting system at the zero-position.**
2. Please transport and store the dispenser only in the original box to avoid damages. Storage temperature: -10° C (+14° F) to +50° C (+122° F).

VIII. Maintenance and Cleaning

1. To guarantee a trouble-free function also for a long term it is necessary to clean the dispenser regularly.
2. Do not leave the dispenser empty without having it cleaned before to avoid blocking of the valve balls! If this happens all the same please put the dispenser for some hours in cold water or clean it by using ultrasound. Once the piston can be moved again, dismantle and remove any residual impurities by cleaning with a normal solution.
3. The dispenser can be left on the bottle all the time if the channel system is filled.
4. When changing the dispensing liquid the dispenser should be cleaned.

IX. Repair Service

1. The POLYFIX® dispenser is mainly maintenance-free and reliable. If there is a problem or a defect and the user himself is not able to solve it, the faulty parts will be exchanged through our dealer or the manufacturer. In case a repair would be necessary, it must be done by the manufacturer or an authorized service station. This does not refer to parts which are available as spare parts, such as piston with finger guard, volume setting system, stopper cap, intake tube, adapters.
2. Before returning the dispenser for repair please make sure that it is perfectly cleaned to avoid any contamination.
3. Please operate the dispenser according to the operating manual to be entitled to the warranty. Repairs must be carried out by the manufacturer or authorized service stations.

X. Technical Datas

1. Volume Range - Subdivisions - Accuracy

Volume range in ml	Graduation in ml	Accuracy A in %	Precision CV in %
2,0 - 10,0	0,2	0,6	0,2
10,0 - 50,0	1,0	0,6	0,2

The datas of accuracy and precision are established according to the German calibration ordinance (PTB) at a room temperature of +20° C (+68° F) by using distilled water. By using chemicals with a different viscosity or density a calibration on an analytical balance is necessary.

2. Packing

POLYFIX® dispensers 10ml and 50ml: single carton, 25,5 x 15,5 x 7 cm.

3. Weight, Dimensions

Volume range in ml	Weight in g appr.	Height in mm
2,0 - 10,0	132	160
10,0 - 50,0	237	220

XI. Equipment

The dispenser comes in a break-proof single carton with 1 PTFE intake tube, 2 PP adapters and 1 operating manual.

XII. Warranty

The POLYFIX® dispenser probably is the most reliable bottle top dispenser in the world. Therefore the manufacturer guarantees for manufacture and original parts for a 12 months period excluding all claims for damage, conversion and other claims. Under this warranty all parts that have been defective due to faulty material or production will be replaced. The manufacturer is not responsible for any damages or problems caused by transport or improper handling. All claims must be reported to the supplier. After a warranty repair the warranty period will not be extended. This warranty is not transferrable and void if repairs are carried out by unauthorized service stations.

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All rights reserved, especially the rights to illustrations and original texts in their whole and in part, for any and all reproduction.